Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Brugsanvisning / Instruksjoner for bruk / Användningsinstruktioner / Käyttöohjeet / Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Inštrukcie na používanie / Használati útmutató / Instruțiuni de folosire / Navodila za uporabo / Οδηγίες χρήσης / Инструкции по применению / Інструкція із застосування / 使用説明 / 使用説明書 / 사용 지침 / تعليمات الاستخدام

# Respireo **SOFT**


## Facial Non Vented

Masque oronasal / Oronasal mask / Mund-Nasen-Maske / Mascariilla oro-Nasal / Maschera oronasale / Máscara oro-nasal / Oronasal maske / Oronasal maske / Oronasal mask / Oronasaalinen kasvomaski / Oronasaal masker / Maska oronasalna / Oronasální maska / Oronazálna maska / Oronasal maszk / Mască oronazală / Oronazalna maska / Στοματορινική μάσκα / Оронозальна маска / Оронозальна маска / 口鼻面罩 / 鼻と口のマスク / 오로나살 마스크 / قناع فموي أنفي


A

**B**


1



2




3


**B**


**D**


1


2


3


**E**

2

3

FRANÇAIS .....	p. 10	ROMÂNĂ.....	p. 50
ENGLISH.....	p. 12	SLOVENŠČINA.....	p. 53
DEUTSCH .....	p. 15	ΕΛΛΗΝΙΚΑ .....	p. 55
ESPAÑOL .....	p. 18	РУССКИЙ.....	p. 58
ITALIANO .....	p. 21	УКРАЇНЕЦЬ.....	p. 61
PORTUGUÊS .....	p. 23	繁體中文.....	p. 64
DANSK.....	p. 26	日本語.....	p. 66
NORSK .....	p. 29	한국어.....	p. 68
SVENSKA.....	p. 31	اللغة العربية.....	p. 70
SUOMI.....	p. 34		
NEDERLANDS.....	p. 37		
POLSKI .....	p. 40		
ČEŠTINA.....	p. 42		
SLOVENČINA.....	p. 45		
MAGYAR.....	p. 48		

## 1 - UTILISATION

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## ENGLISH

## 1 - INTENDED USE

Respireo SOFT Facial Non Vented is a reusable oronasal mask without calibrated exhalation holes intended to be used at home or in a hospital setting by multiple adult patients (weighing more than 30 kg) under Non invasive Positive-Pressure Ventilation (NPPV).

 **WARNINGS**

- The mask should only be used with a ventilator fitted with an active exhalation valve and/or with a respiratory circuit with an exhalation valve recommended by your doctor or respiratory therapist.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.

- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, or exceeding the recommended number of regeneration cycles, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 2 - MASK COMPONENTS (section A)

a) cushion; b) frame; c) blue elbow connector; d) headgear; e) buckle.

## 3 - ASSEMBLING AND WEARING THE MASK (section B)

- Holding the mask gently in close contact to the face, wear the headgear (d) over the head making sure the cushion (a) adheres to the oronasal profile (B.1).
- Attach the two lower straps of the headgear (d) to the hooks of the frame (b) so that they pass under the ears (B.2).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.3.
- Connect the blue elbow connector (c) to the equipment tube (B.4), then turn on the equipment at the prescribed pressure.
- Lie down and, if necessary, adjust the headgear again (B.5).

## 4 - USING THE MASK

To properly remove, disassemble and reassemble the mask, refer to the respective illustrations in section C - D - E.

## 5 - CLEANING THE MASK

For daily cleaning, the mask must be disassembled in its components, each of which (except for the headgear) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

**To clean the headgear, perform the same procedure weekly.**

For a deeper decontamination of the mask, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website [medicaldevice.airliquide.com](http://medicaldevice.airliquide.com) or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

## PRECAUTIONS

- Remove the mask after disconnecting at least one of the lower headgear straps.
- Do not open the straps to avoid having to adjust the headgear again the next time you use it.
- Do not disassemble the blue elbow connector from the frame under any circumstances, even when cleaning the mask.
- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

## 6 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS

The mask can be regenerated for up to 20 cycles by one of the following processes:

- high level thermal disinfection;
- high level chemical disinfection;
- autoclave sterilization at 121°C for a maximum of 30 minutes.

For full detailed mask regeneration instructions, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website [medicaldevice.airliquide.com](http://medicaldevice.airliquide.com) or at your dealer.

### PRECAUTIONS

- *The mask must always be regenerated before use on a new patient.*
- *The headgear must be replaced.*

## 7 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> <li>- Flow generator turned off or not working.</li> <li>- Breathing circuit ventilator not inserted correctly.</li> </ul>	<ul style="list-style-type: none"> <li>- Switch on or adjust the flow generator.</li> <li>- Connect the breathing circuit ventilator to the flow generator and the mask.</li> </ul>
Noisy and annoying air leaks.	<ul style="list-style-type: none"> <li>- The mask has not been assembled correctly.</li> <li>- The mask is not positioned correctly.</li> <li>- The size of the mask may not be correct.</li> </ul>	<ul style="list-style-type: none"> <li>- Disassemble the mask and then reassemble it as described in the user manual.</li> <li>- Reposition the mask on your face and adjust the headgear straps as described in the user manual.</li> <li>- Contact your doctor or therapy technician to find the right size.</li> </ul>

## 8 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polycarbonate
Frame (b)	Polycarbonate
Blue elbow connector (c)	Polycarbonate
Headgear (d)	Nylon/UBL/Polyurethane
Buckle (e)	Nylon

## 9 - TECHNICAL DATA

<b>Regulatory information</b>	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 5356-1
<b>Operating pressure range</b>	4 cmH <sub>2</sub> O / 32 cmH <sub>2</sub> O
<b>Operating temperature range</b>	+5 °C / +40 °C
<b>Dead space</b>	Size S: 281 mL Size M: 293 mL Size L: 313 mL
<b>Dimensions (without headgear)</b>	Size S: 190 x 110 x 110 mm Size M: 190 x 110 x 110 mm Size L: 190 x 110 x 110 mm
<b>Weight</b>	Size S: 117 g Size M: 118 g Size L: 120 g

## 10 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case

of product defects can be applied by contacting your dealer.  
For more information on the rights under the warranty, please contact your retailer.

**COPYRIGHT**

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

**TECHNICAL CHARACTERISTICS UPDATES**

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market.

If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.