FRANÇAIS

1. UTILISATION

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ENGLISH

1. INTENDED USE

Respiro Hospital F2 Non Vented Disposable with anti-asphyxia valve is an oronasal mask without functional exhalation orifices and with an anti-suffocation valve, intended to be used for short time periods (less than 7 days) in a hospital setting by a single adult patient (weighing over 30kg) in controlled non-invasive positive pressure ventilation (NIPPV).

⚠️ WARNINGS

- The mask should only be used with a ventilator fitted with an active exhalation valve and/or with a respiratory circuit with an exhalation valve recommended by your doctor or respiratory therapist.
- The anti-asphyxia valve fitting ensures patient safety even in case of failure of the positive pressure equipment. Do not use the mask if the valve is damaged or distorted.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- The device is a single use device: the reuse on multiple patients may result in transmission of infections, serious harm or even the death of the patient.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2. TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask	- Flow generator turned off or not working - Breathing circuit ventilator to the flow generator not inserted correctly	- Switch on or adjust the flow generator - Connect the breathing circuit ventilator to the flow generator and the mask
Noisy and annoying air leaks	- The mask has not been assembled correctly - The mask is not positioned correctly	- Disassemble the mask and reassemble it as described in the user manual - Reposition the mask on your face and adjust the headgear straps as described in the user manual - Contact your doctor or therapy technician to find the right size
- The size of the mask may not be correct		

3. DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone
Frame (b)	Polycarbonate
Forehead support (c)	Silicone
Anti-asphyxia valve fitting with yellow connector (d)	Polycarbonate + Silicone
Headgear (e)	Nylon/UBL/ Polyuréthane
Clip (f)	Nylon

4. TECHNICAL DATA

Regulatory information	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 5356-1
------------------------	--

Operating pressure range	4 cmH ₂ O / 32 cmH ₂ O
Operating temperature range	+5 °C / +40 °C
Opening pressure of the anti-asphyxia valve	0,55 cmH ₂ O
Closing pressure of the anti-asphyxia valve	0,95 cmH ₂ O
Dead space	S: 287 ml M: 297 ml L: 328 ml
Dimensions (without headgear)	S: 200x100x115 mm M: 200x105x115 mm L: 210x100x115 mm
Weight	S: 138 g M: 141 g L: 142 g

5. WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components. The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use. Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

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